

## Noncovalent Protein Internalization into Mammalian Cells

### Introduction

Introduction of active and properly functioning proteins into mammalian cells has proved a useful tool in the treatment of various diseases. It is now known that proteins can be delivered into both mammalian and plant cells using arginine-rich intracellular delivery (AID) peptides in a noncovalent manner. This technique allows for the quick and simple internalization of functioning proteins into cells. This protocol outlines the use of arginine-rich intracellular delivery (AID) peptides to deliver fluorescent proteins or  $\beta$ -galactosidase enzymes into animal or plant cells, providing a useful strategy to introduce active proteins in cells and tissues *in vivo*. In this protocol, X-Gal is used to provide a visual indicator since it can be hydrolyzed by  $\beta$ -galactosidase, yielding 5,5'-dibromo-4,4'-dichloro-indigo-2, a blue product.

### Materials


- $\beta$ -Gal protein
- AID peptides in PBS
- MCF7 cells
- PBS (Phosphate-Buffered Saline) (GoldBio Catalog # [P-271](#))
- X-Gal (GoldBio Catalog # [X4281C](#))

### Method

In this procedure, cells are treated with  $\beta$ -Gal (control) or  $\beta$ -Gal/AID (experimental), followed by X-Gal treatment. After washing the X-gal away, cells can then be observed under a microscope.

1. For the control sample, aspirate medium and treat cells with 0.5 $\mu$ M  $\beta$ -Gal (5  $\mu$ g  $\beta$ -Gal in a sample volume of 70.8  $\mu$ l) for 15 minutes at room temperature.
2. For experimental samples, prepare a solution containing 0.5 $\mu$ M  $\beta$ -Gal and 12.1 $\mu$ M AID peptides (1:24 ratio), vortex for 10 seconds and incubate at room temperature for 20 minutes. Treat experimental cells with this  $\beta$ -Gal/AID solution at room temperature for 15 minutes.
3. Aspirate the  $\beta$ -Gal or  $\beta$ -Gal/AID solution from the control or experimental samples, respectively, and wash with 1 ml PBS three times.
4. Treat cells with 2 mg X-Gal in 1 ml PBS and incubate at 37°C for 8 hours.

5. Aspirate the X-Gal solution and wash cells with 1 ml PBS three times.
6. Observe cells under microscope. Determine internalization by localizing blue product in cells as shown below.



**Fig. 2. Confocal microscopy of noncovalent protein internalization via AID peptides in animal cells.** (E) Image of cells incubated with b-Gal in bright field. Cells were treated with 0.5 $\mu$ M b-Gal followed by X-gal treatment. Image was shown at a magnification of 1000x. (F) Image of cells treated with R9-plus-b-Gal mixtures in bright field. Cells were treated with 0.5 $\mu$ M b-Gal pre-mixed with 12.1 $\mu$ M of the R9 peptide followed by X-gal treatment.

### Tips

- MCF7 cells are maintained in RPMI 1640 medium supplemented with 10% heat-inactivated bovine serum and 1X penicillin/streptomycin/amphotericin B.

### Associated Products

- [PBS \(GoldBio Catalog # P-271\)](#)
- [X-Gal \(GoldBio Catalog # X4281C\)](#)

### References

Wang, Y. H., Chen, C. P., Chan, M. H., Chang, M., Hou, Y. W., Chen, H. H., Hsu, H.R., Liu, K., and Lee, H. J. (2006). Arginine-rich intracellular delivery peptides noncovalently transport protein into living cells. *Biochemical and biophysical research communications*, 346(3), 758-767.