

Nucleic Acids Extraction - Nuclease Inactivation Utilizing Proteinase K

Introduction

Proteinase K is used in the extraction of nucleic acids. After cell lysis, nucleases are released that degrade DNA and RNA. Proteinase K can effectively inactivate these nucleases by digesting them. Proteinase K can also help by digesting other proteins that are capable of contaminating your nucleic acid sample. Proteinase K can be used in most any DNA or RNA extraction protocol after cell lysis, but before extraction. The use of proteinase K is described below in a general nucleic acid extraction protocol.

Materials

- Proteinase K ([Proteinase K, GoldBio Catalog # P-480](#) [CAS 39450-01-6, mw.=28.9 kDa])
- Any other materials that are used in your extraction protocol

Method

Cell Lysis

1. Lyse cells using your lysis method of choice.

Note: Proteinase K is compatible with [guanidinium chloride](#), [guanidinium thiocyanate](#), [urea](#), [iodoacetate](#), [citrate](#), [sodium dodecyl sulfate \(SDS\)](#), [Triton X-100](#), [Tween 20](#) and [EDTA](#).

Nuclease Digestion

2. Add Proteinase K to the lysate (See [Proteinase K Stock Solution](#) protocol).

Note: The final concentration of Proteinase K in solution should be 50-400 µg/ml.

3. Incubate sample at 55°C for 1- 3 hours.

Note: Samples may be digested overnight to be sure of complete nuclease digestion.

Note: Proteinase K can then be inactivated by heating to 95°C for 10 minutes after digestion.

Extraction

4. Extract the desired nucleic acid using your extraction method of choice.

References

Hilz, H., Wieggers, U., & Adamietz, P. (1975). Stimulation of proteinase K action by denaturing agents: application to the isolation of nucleic acids and the degradation of 'masked' proteins. *European Journal of Biochemistry*, 56(1), 103-108.

Goldenberger, D., Perschil, I., Ritzler, M., & Altwegg, M. (1995). A simple "universal" DNA extraction procedure using SDS and proteinase K is compatible with direct PCR amplification. *Genome Research*, 4(6), 368-370.

McGookin, R. (1984). RNA extraction by the proteinase K method. In *Nucleic Acids* (pp. 109-112). Humana Press.

Maloy, S. R. (1990). *Experimental techniques in bacterial genetics*. Jones & Bartlett Learning.