

Stock Solution



TD-S Revision 2.0

Creation Date: 8/1/2014
Revision Date: 4/15/2019

Picloram 10 mg/ml Stock Solution

Instructions

1. Add 1 g of Picloram [[GoldBio Catalog # P0914](#), (CAS 1918-02-1, MW = 241.5 g/mol)] to a 100 ml volumetric flask, beaker or other glass container.
2. Add 75 ml of Acetonitrile or DMSO to dissolve the powder. Stirring may be required to completely dissolve the Picloram. Once completely dissolved, bring to 100 ml final volume with Acetonitrile or DMSO.
3. The stock solution may be stored at -20 C.