



Product Specification

Product Name: Murine LIF
Cat #: 1330-07
CAS #:
Storage: Store desiccated at -20°C.

Test	Specifications
Characteristics	White Freeze-dried Powder
SDS-PAGE Electrophoresis Purity	>98.0%
HPLC Purity	>98.0%
Specific Biological Activity	>1.0 × 10 ⁸ IU/mg M1 myeloid leukemic cells induction = <0.01 ng/ml
Molecular Weight	19.9±2.0 KDa
Bacterial Endotoxins	<1.0 EU/μg rMuLIF

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com