



Product Specification

Product Name: Recombinant HBV Surface Antigen-preS1
Cat #: 6820-01
CAS #:
Storage: Store desiccated at -20°C.

Test	Specifications
Characteristics	White Freeze-dried Powder
SDS-PAGE Electrophoresis Purity	>95%
HPLC Purity	>95%
Molecular Weight	12.6±1.3 kDa (dimer)
Bacterial Endotoxins	<1.0 EU/μg rHBsAg-preS1

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com