



Product Specification

Product Name: 4% Agarose Beads, Standard

Cat #: A-182

CAS #:

Storage: Store at room temperature.

Test	Specifications
Bead Geometry	Spherical
Bead Size Range	50 -150 μ m
Selectivity (K_{av}) - $M_p=1.000.000$ - $M_p=100.000$ - $M_p=20.000$	0.15 - 0.39 0.50 - 0.66 0.75 - 0.87
Column Flow Rate Test ($p/L=15$; cm/h)	70-140

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com