



Product Specification

Product Name: 4% Agarose Beads, Standard, Crosslinked
Cat #: A-192
CAS #: N/A
Storage: Store at room temperature.

Test	Specifications
Bead Geometry	Spherical
Bead Size Range	50 -150 µm
pH Stability	Passes
Chemical Stability	Passes
Cytochrome C Binding	Passes
Autoclavable (20 min. at 120°C)	Passes
Column Flow Rate Test	Passes

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com