



Product Specification

Product Name: 6% Agarose Beads, Standard, Crosslinked
Cat #: A-193
CAS #:
Storage: Store at room temperature.

Test	Specifications
Bead Geometry	Spherical
Bead Size Range	50-150 μ m
Flow Rate (p/L=15)	100-200 cm/h
Microbial Contamination (Colony forming units / ml suspension)	$\leq$ 100 units

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com