



Product Specification

Product Name: L-Arabinose
Cat #: A-300
CAS #: 87-72-9
Storage: Store at room temperature.

Test	Specifications
Appearance	White or almost white crystalline powder
Coliform	≤3 MPN/g
Specific rotation	+102° to +105°
Melting point	154-160°C
Loss on Drying	≤0.5%
Ash	≤0.5%
Chloride	≤50 mg/kg
Sulfate	≤50 mg/kg
Heavy metals (Pb)	≤0.5 µg/g
Arsenic (As)	≤1.0 µg/g
Total plate count	≤1000 cfu/g
Assay	99.0 - 102.0%

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com