



Product Specification

Product Name: D-Arabitol
Cat #: A-310
CAS #: 488-82-4
Storage: Store at room temperature.

Test	Specifications
Color	white to off-white
Form	crystalline powder
Color of Solution (10% water)	Clear or colorless
Water Content	max. 1.0%
Purity	min. 99.0 area-%
Specific rotation (c = 5 in borax solution 8%)	10.0 to 12.0°
Melting Point	100.0-104.0°C
FT-IR Spectrum	corresponds to reference spectrum

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com