



Product Specification

Product Name: L-Arabitol
Cat #: A-320
CAS #: 7643-75-6
Storage: Store desiccated at 4°C. Protect from light.

Test	Specifications
Color	white to off-white
Form	crystalline powder
Color of Solution (10% in water)	colorless
Turbidity (1% in water)	max. 4.0 NTU
Water Content (K. Fischer)	max. 1.0%
Melting Point	100 to 104°C
Specific Optical Rotation [alpha] _D (c = 5 in borax solution 8%, 20°C)	-13.0 to -11.0
Purity (HPLC)	min. 98.0 area-%
FT-IR Spectrum	corresponds to reference spectrum

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com