



Product Specification

Product Name: Arachidoyl-DL-carnitine chloride
Cat #: A-327
CAS #: 149116-07-4
Storage: Store at -20°C. Protect from light.

Test	Specifications
Color	white to off-white
Form	powder
Color of Solution (1% in methanol)	clear, colorless
Water Content (K. Fischer)	max. 1.0%
TLC (silica gel, 50 mg in 2 ml methanol; 1-Butanol/acetic acid/water 2:1:1; iodide)	corresponds to reference material
Assay (titration with 0.1 mol/L Silver Nitrate)	min. 95.0%
FT-IR Spectrum	corresponds to reference spectrum

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com