



## Product Specification

**Product Name:** Bovine Serum Albumin (BSA), Cohn Analog  
**Cat #:** A-422  
**CAS #:** 9048-46-8  
**Storage:** Store at 4°C.

| Test                          | Specifications                                               |
|-------------------------------|--------------------------------------------------------------|
| Appearance                    | White to yellow with tan to green cast<br>lyophilized powder |
| Solubility<br>(4% in water)   | Clear to slightly hazy                                       |
| Protein (db)                  | ≥96.0%                                                       |
| Purity                        | ≥96.0%                                                       |
| Moisture                      | ≤5.0%                                                        |
| pH<br>(10% solution in water) | 6.5 to 7.5                                                   |
| Ash<br>(Residue on Ignition)  | <3.0%                                                        |
| Endotoxin                     | as reported                                                  |

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

**Gold Biotechnology, Inc.**

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: [contactgoldbio86@goldbio.com](mailto:contactgoldbio86@goldbio.com)