



Product Specification

Product Name: Bovine Serum Albumin (BSA), Sulfhydryl Blocked (NZ)
Cat #: A-424
CAS #: 9048-46-8
Storage: Store at 4°C.

Test	Specifications
Appearance	White to yellow with tan to green cast lyophilized powder
Protein (db)	≥98.0%
Purity	≥98.0%
Moisture	≤5.0%
pH (7% in Water)	6.5 - 7.5
Sodium	≤15.0 mg/g
BIgG	None Detected
Chloride	≤6.0 mg/g
Protease	<0.001 Units/mg
Free Sulfhydryl Content	≤0.1 mol/mol Albumin
Absorbance	≤0.100
Asepsis	None Detected
BVD Virus	None Detected
Mycoplasma	None Detected
Pathogen Inducing Cytopathy	None Detected
Pathogen Inducing Haemadsorption	None Detected

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com