



Product Specification

Product Name: 8-Acetoxypyrene-1,3,6-trisulfonic acid trisodium salt
Cat #: A-850
CAS #: 115787-83-2
Storage: Store at -20°C.

Test	Specifications
UV Wavelength (0.02 g/L in 90 mmol/L tris buffer pH 8, lambda max.)	367-373 nm
Molar Extinction Coefficient (0.02 g/L in 90 mmol/L tris buffer pH 8, at lambda max.)	≥25000 L/(mol cm)
Color	slightly yellow to yellow
Form	crystalline powder
Color of Solution (1% in water)	colorless to yellow to green
Turbidity (1% in water)	≤4.0 NTU
Purity (HPLC)	≥98.0 area-%

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com