



Product Specification

Product Name: Bis-Tris
Cat #: B-020
CAS #: 6976-37-0
Storage: Store at room temperature.

Test	Specifications
Appearance	White to off-white crystalline powder
Assay	99.0-101.0%
Water Content (K. Fisher)	≤1.0%
pH	9.5-11.0
Solubility (0.1M in aqueous solution)	Clear, colorless solution
UV absorbance - 280 nm (0.1M in aqueous solution)	≤0.040
UV absorbance - 340 nm (0.1M in aqueous solution)	≤0.020

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com