



Product Specification

Product Name: Lapis-3-Acetate
Cat #: B-650
CAS #: 17357-14-1
Storage: Store at -20°C. Protect from light.

Test	Specifications
Color	white to off-white
Form	crystalline powder
Color of Solution (5% in ethanol)	colorless to almost colorless
Turbidity (0.5% in ethanol)	max. 4.0 NTU
Water Content (K. Fischer)	max. 0.5%
Melting Point	130 - 134°C
Purity (HPLC)	min. 98.0 area-%
FT-IR Spectrum	corresponds to reference spectrum

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com