



Product Specification

Product Name: Lapis-Caprylate
Cat #: B-660
CAS #: 133950-69-3
Storage: Store at -20°C. Protect from light.

Test	Specifications
Color	white to off-white
Form	crystalline powder
Color of Solution (5% in acetone)	colorless
Turbidity (0.5% in acetone)	max. 2.0 NTU
Water Content (K. Fischer)	max. 1.0%
Purity (TLC) (30 mg/mL in acetone; toluene/acetone 5:1; UV 254 nm)	min. 99 area-%
TLC (30 mg/mL in acetone; toluene/acetone 5:1; UV 254 nm)	corresponds to reference material
FT-IR Spectrum	corresponds to reference spectrum

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com