



Product Specification

Product Name: Bis-Acrylamide
Cat #: B-900
CAS #: 110-26-9
Storage: Store at 4°C.

Test	Specifications
Appearance	White crystalline powder
Conductivity (2% in water)	≤5 μS/cm
pH (1% in 0.1M NaCl)	≥5.0
Absorbance (290 nm, 1% in water)	≤0.2
Acrylic Acid	≤0.001%
DNase	None Detected
Rnase	None Detected
Protease	None Detected
Assay	≥99.0%
Heavy Metals (Pb)	≤1.5 ppm
Water Content	≤1.0%

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com