



Product Specification

Product Name: Bis-Acrylamide, 2% Solution
Cat #: B-901
CAS #: 110-26-9
Storage: Store at 4°C.

Test	Specifications
Appearance	Clear and colorless solution
Conductivity (2% in Water)	≤5 μS/cm
Acrylic Acid	≤0.001%
Assay (powder, HPLC)	≥99.0%

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com