



Product Specification

Product Name: Caffeic acid phenethyl ester
Cat #: C-070
CAS #: 104594-70-9
Storage: Store desiccated at -20°C. Protect from light.

Test	Specifications
Color	white to yellow
Form	crystalline powder
Color of Solution (5% in acetone)	colorless to yellow
Turbidity (0.5% in acetone)	max. 2.0 NTU
Water Content (K. Fischer)	max. 0.5%
Purity (HPLC)	min. 99.0 area-%
Impurity (phenethyl alcohol; HPLC)	max. 0.50%
FT-IR Spectrum	corresponds to reference spectrum

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com