



Product Specification

Product Name: Carbenicillin (Disodium) 100 mg/mL Solution
Cat #: C-103-SL
CAS #: 4800-94-6
Storage: Store at -20°C.

Test	Specifications
Appearance	Clear solution
Water	≤6.0%
pH	5.5-7.5
Iodine absorbing substances	≤8.0%
Penicillin G (Benzylpenicillin sodium)	≤5.0%
Assay (anhydrous basis)	≥88.0%
Potency (on Dried Basis)	≥770 µg/mg
Solution	100 mg/ml
Sterile Filtered	Complies

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com