



Product Specification

Product Name: Dehydroascorbic acid
Cat #: D-080
CAS #: 490-83-5
Storage: Store desiccated at -20°C. Protect from light.

Test	Specifications
Color	white to light yellow or light brown
Form	powder
Color of Solution (1% in water)	colorless to slightly yellow or slightly brown
Turbidity (1% in water)	≤4.0 NTU
Water Content (K. Fischer)	≤1.0%
Sodium Hydroxide Titration (0.1 mol/L)	97.0 to 103.0%
FT-IR Spectrum	corresponds to reference spectrum

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com