



Product Specification

Product Name: DAP-HCl
Cat #: D-130
CAS #: 54897-59-5
Storage: Store at 4°C. Protect from light.

Test	Specifications
Color	white to off-white
Form	crystalline powder
Color of Solution (5% in water)	colorless
Turbidity (0.5% in water)	max. 2.0 NTU
Water Content (K. Fischer)	max. 0.5%
Melting Point	223 - 227°C
Argentometric Titration	99.0 - 101.0%
TLC (50 mg in 1.5 ml water; 2-propanol/ammonia 2:1; ninhydrin)	corresponds to reference material
FT-IR Spectrum	corresponds to reference spectrum

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com