



Product Specification

Product Name: DTNB (Ellman's Reagent)
Cat #: D-250
CAS #: 69-78-3
Storage: Store at 4°C. Protect from light and heat.

Test	Specifications
Color	white to slightly yellow to yellow
Form	powder
Water Content (K. Fischer)	≤1.0%
Purity (HPLC)	≥98.0 area-%
¹ H-NMR Spectrum	corresponds to assigned structure
¹³ C-NMR Spectrum	corresponds to assigned structure

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com