



Product Specification

Product Name: DNase I, Bovine Pancreas >2000 Kunitz U/mg
Cat #: D-301
CAS #: 9003-98-9
Storage: Store desiccated at -20°C.

Test	Specifications
Appearance	White or almost white crystalline powder
Potency	≥2000 U/mg

One Kunitz Unit will produce an A260 of 0.001 per minute at pH 5.0 at 25°C using DNA, type I or III as the substrate.

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com