



Product Specification

Product Name: DNase I, Bovine Pancreas (Lysate Tested), >300 Kunitz U/mg
Cat #: D-303
CAS #: 9003-98-9
Storage: Store desiccated at -20°C.

Test	Specifications
Appearance	White or almost white powder
Enzyme activity (on Dried Basis)	≥300 Kunitz u/mg
Loss on Drying	≤5.0%

Unit Definition: One Kunitz unit will produce A260 of 0.001 per minute at pH 5.0 at 25°C using DNA, type I or III as the substrate.

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com