



Product Specification

Product Name: dNTP Mix
Cat #: D-900
CAS #:
Storage: Store at -20°C.

Test	Specifications
Appearance	Clear, aqueous solution
Purity (HPLC)	≥99.0%
Performance - Taq DNA Polymerase - Pfu DNA Polymerase - Q-PCR - Long PCR (>30 kb)	Passes Passes Passes Passes
DNase	None
RNase	None
Phosphatase	None

Avoid freeze/thaw cycles or exposure to frequent temperature changes.

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com