



Product Specification

Product Name: Sodium D-glucuronate
Cat #: G-140
CAS #: 207300-70-7
Storage: Store at room temperature.

Test	Specifications
Color	white to off-white
Form	powder
Solubility (50 mg/mL in H ₂ O)	clear colorless to yellow solution
Water Content (K. Fischer)	≤8.5%
Specific Optical Rotation [alpha] _D (c = 2 in water)	20.0 to 23.0
Titration	≥98.0%
FT-IR Spectrum	corresponds to reference spectrum

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com