



Product Specification

Product Name: Guanidine thiocyanate

Cat #: G-210

CAS #: 593-84-0

Storage: Store desiccated at room temperature. Protect from light.

Test	Specifications
Appearance	white crystalline powder
Solubility	Clear, colorless solution
pH (1M aqueous solution)	4.5-7.0
Water Content	≤0.5%
Melting Point	115.0-120.0°C
Assay	≥99.0%
UV-Absorption (6M water at 410 nm)	≤0.050
UV-Absorption (6M water at 300 nm)	≤0.100
Heavy Metals	<5 ppm

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com