



Product Specification

Product Name: Nickel Agarose Beads (High Density)
Cat #: H-320
CAS #:
Storage: Store at 4°C.

Test	Specifications
Bead Geometry	Spherical
Bead Size	Standard - 50-150 μm
Matrix	6% crosslinked agarose bead
Binding Capacity	20-40 $\mu\text{mol Ni}^{2+}$ / ml gel

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com