



Product Specification

Product Name: N-gamma-Hydroxy-L-arginine, acetate salt
Cat #: H-360
CAS #: 53598-01-9
Storage: Store desiccated at -20°C.

Test	Specifications
Color	off-white to beige
Form	lyophilisate
Specific Optical Rotation [alpha] _D (c = 0.5 in water; 20°C)	8.0 - 12.0
Purity (TLC) (10 mg/mL in water; butanol/acetic acid/water 3:1:1, ninhydrin)	≥95.0 area-%
FT-IR Spectrum	corresponds to reference spectrum

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com