



Product Specification

Product Name: Illumination™ Firefly & Renilla Luciferase Enhanced Assay Kit
Cat #: I-920
CAS #:
Storage: Store the kit at -80°C.

Test	Specifications
Luciferase Reporter Assay (with purified recombinant firefly luciferase)	Titration curve meets expected standard range
Luciferase Reporter Assay (with purified recombinant <i>Renilla</i> luciferase)	Titration curve meets expected standard range
Luciferase Reporter Assay (with cell lysate transfected with firefly luciferase vector)	Titration curve meets expected standard range
Luciferase Reporter Assay (with cell lysate transfected with <i>Renilla</i> luciferase vector)	Titration curve meets expected standard range

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com