



Product Specification

Product Name: L-Valine
Cat #: L-090
CAS #: 72-18-4
Storage: Store at room temperature.

Test	Specifications
Color	white
Form	powder
Solubility (50 mg/mL in water)	colorless, clear
Specific Optical Rotation $[\alpha]_D$ (c = 8 in 6 mol/L HCl at 25°C)	26.6 to 28.4
Perchloric Acid Titration (0.1 mol/L)	≥98%
FT-IR Spectrum	corresponds to reference spectrum

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com