



Product Specification

Product Name: D-Luciferin, Sodium Salt
Cat #: LUCNA
CAS #: 103404-75-7
Storage: Store desiccated at -20°C. Protect from light.

Test	Specifications
Color	slightly yellow to yellow
Form	powder
Color of Solution (1% in water)	slightly yellow to yellow
Turbidity (1% in water)	≤4.0 NTU
Water Content (K. Fischer)	≤8.0%
TLC (50 mg in 5 ml water; ethanol/ammonia 20:1; UV 254 nm)	Corresponds to reference material
Purity (HPLC)	≥99.0 area-%
UV Wavelength (c=0.01 mg/mL in 0.5 M Na ₂ CO ₃)	383-387 nm
Molar Extinction Coefficient	≥17500.0 L/(mol cm)
UV Wavelength (c=0.01 mg/mL in 0.5 M Na ₂ CO ₃)	282-286 nm
Molar Extinction Coefficient	≥6000.0 L/(mol cm)
FT-IR Spectrum	Corresponds to reference spectrum

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com