



Product Specification

Product Name: Phenyl- β -D-glucuronide
Cat #: P-140
CAS #: 17685-05-1
Storage: Store at -20°C. Protect from light.

Test	Specifications
Color	white to off-white
Form	crystalline powder
Water Content (K. Fischer)	$\leq 8.0\%$
Purity (TLC)	$\geq 98.0\%$
Specific Optical Rotation $[\alpha]_D$ (c =1 in water)	-92.0 to -82.0°
FT-IR Spectrum	corresponds to reference spectrum

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com