



Product Specification

Product Name: D-Proline
Cat #: P-220
CAS #: 344-25-2
Storage: Store desiccated at room temperature.

Test	Specifications
Appearance	White to off-white powder
Purity (Titration)	≥98.0%
Loss on drying	≤0.5%
Specific Optical Rotation [α] _D (c = 4 in water, 20°C)	83.0 to 86.5

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com