



Product Specification

Product Name: DL-Proline
Cat #: P-230
CAS #: 609-36-9
Storage: Store desiccated at room temperature.

Test	Specifications
Appearance	White to off-white powder
Purity	≥98.5%
Loss on drying	≤0.5%
Residue on ignition	≤0.4%
Heavy metals	≤15 ppm
Chloride	≤0.2%
Iron	≤30 ppm
Sulfate	≤0.2%

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com