



Product Specification

Product Name: Polyethylene Glycol (PEG) 3350
Cat #: P-290
CAS #: 25322-68-3
Storage: Store desiccated at room temperature.

Test	Specifications
Color	White to off-white granular powder
Form	powder or crystalline powder
pH (5%, Water at 25°C)	4.5-7.5
FT-IR Spectrum	corresponds to reference spectrum

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com