



Product Specification

Product Name: Protein G HTC Agarose Beads

Cat #: P-430

CAS #:

Storage: Store at 4°C.

Test	Specifications
Bead Geometry	Spherical
Bead Size	Standard - 50-150 μ m
Matrix	4% crosslinked agarose bead
Binding Capacity	10 - 30 mg IgG _h / ml gel

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com