



Product Specification

Product Name: Protein G Magnetic Agarose Beads
Cat #: P-434
CAS #:
Storage: Store at 4°C. Do not freeze.

Test	Specifications
Bead Geometry	Spherical
Bead Size	Standard - 50-100 µm
Matrix	4% Crosslinked Agarose
Static Binding Capacity	~20 mg human IgG/ml gel

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com