



Product Specification

Product Name: Puromycin Dihydrochloride
Cat #: P-600
CAS #: 58-58-2
Storage: Store desiccated at -20°C.

Test	Specifications
Appearance	White to off-white powder
Solubility (50 mg/ml in water)	Clear, colorless to yellowish solution
Melting Point	160-182°C
Purity (HPLC)	≥98.0%
TLC (NP, CH ₂ Cl ₂ , Methanol 7:3)	≥99.0%
TLC (RP C18, methanol, H ₂ O 9:1)	≥99.0%
Water Content	≤12.0%

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com