



Product Specification

Product Name: 6% HTC Agarose Beads, Standard
Cat #: R-201
CAS #:
Storage: Store at room temperature.

Test	Specifications
Bead Geometry	Spherical
Bead Size	50-150 μm
Matrix	6% agarose beads
Cationic Interaction (Cytochrome C / 1 ml gel bed)	<8 μg
Column Flow Rate	>6.8 ml/min
Maximum Column Flow Rate/Pressure	<65 ml/min

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com



Product Specification

Product Name: 6% HTC Agarose Beads, Standard
Cat #: R-201
CAS #:
Storage: Store at room temperature.

Test	Specifications
Bead Geometry	Spherical
Bead Size	50-150 μm
Matrix	6% agarose beads
Cationic Interaction (Cytochrome C / 1 ml gel bed)	<8 μg
Column Flow Rate	>6.8 ml/min
Maximum Column Flow Rate/Pressure	<65 ml/min

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com