



Product Specification

Product Name: Streptomycin Sulfate 100 mg/mL Solution
Cat #: S-150-SL
CAS #: 3810-74-0
Storage: Store at -20°C.

Test	Specifications
Appearance	A white or almost white powder
Methanol (GC)	≤0.3%
Streptomycin B	≤3.0%
Loss on drying	≤7.0%
Sulfated ash	≤1.0%
Sulfate	18.0 to 21.5%
Colormetric test	≥90.0%
pH	4.5 - 7.0
Assay (dried substance)	≥720 U/mg
Solution	100 mg/ml
Sterile Filtered	Complies

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com