



Product Specification

Product Name: Trypsin 1:250
Cat #: T-161
CAS #: 9002-07-7
Storage: Store at -20°C.

Test	Specifications
Appearance	White to off-white powder
Trypsin Activity	≥250 U/mg
Loss on Drying	≤5.0%
Solubility	Soluble, lightly opalescent
Microbial Limits - Staphylococcus aureus - Pseudomonas aeruginosa - Salmonella species	Negative Negative Negative

Conforms to USP Grade Specifications

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com