



Product Specification

Product Name: Trypsin Inhibitor, Soybean, Purified >10,000 U
Cat #: T-166
CAS #: 9035-81-8
Storage: Store desiccated at -20°C.

Test	Specifications
Appearance	Off-white to beige lyophilized powder
Activity (material) (BAEE units of inhibition)	≥10,000 U/mg
Solubility	Complies
Protein Concentration	0.872 mg protein/material

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com