



Product Specification

Product Name: TES
Cat #: T-785
CAS #: 7365-44-8
Storage: Store at room temperature.

Test	Specifications
Appearance	White crystalline powder
Assay (Titration)	99.0 - 101.0%
Water Content (K. Fischer)	≤0.5%
Solubility (1M aqueous)	Clear, colorless solution
pH (1M aqueous)	4.0-5.0
Melting Point	220 - 225°C
Absorbance (A260)	≤0.050
Absorbance (A280)	≤0.030
Heavy Metals	≤5 ppm
Loss on drying	≤1.0%

GOLD BIOTECHNOLOGY, INC. provides this Product Specification in good faith. Specifications are subject to change.

Gold Biotechnology, Inc.

1328 Ashby Rd.

St. Louis, MO 63132

P: (314)890-8778

F: (314)890-0503

E: contactgoldbio86@goldbio.com