

**AFFINITY HIS-TAG PURIFICATION
 TECHNICAL SPECIFICATIONS
 Aminoethyl Agarose BULK resins**

BEAD (Geometry, size)	Spherical, Standard: ~ 50-150 µm		
CROSS-LINKED	Yes		
AUTOCLAVABLE	121°, 30'		
MATRIX ACTIVE GROUPS	Amino groups		
DENSITY	VERY LOW Density	LOW Density	HIGH Density
AGAROSE %	4%	6%	4% 6%
ACTIVATION DEGREE (µm diaminoethyl / mL gel)	3-6	15-25	40-60 40-60
CATALOG #	A-500	A-501	A-499 A-503
ANTIMICROBIAL AGENT	0.5 M NaCl containing 0.02% thimerosal or 20% Ethanol		
STORAGE TEMPERATURE	2 - 8°C		

For laboratory use only. Not for use in diagnostic or therapeutic procedures.

Gold Biotechnology

St. Louis, MO

Web: www.goldbio.com

Ph: (314) 890-8778

email: contactgoldbio86@goldbio.com