

**AFFINITY HIS-TAG PURIFICATION
 TECHNICAL SPECIFICATIONS
 Glyoxal Agarose BULK resins**

BEAD (Geometry, size)	Spherical, Standard: ~ 50-150 µm			
CROSS-LINKED	Yes			
AUTOCLAVABLE	121°, 30'			
MATRIX ACTIVE GROUPS	Agarose with some diols oxidized to aldehydes			
DENSITY	LOW Density		HIGH Density	
AGAROSE %	4%	6%	4%	6%
ACTIVATION DEGREE (µm glyoxal / mL gel)	15-25	15-25	40-60	40-60
BINDING CAPACITY⁽¹⁾ (mg BSA / mL gel)	5-10	5-10	15-20	20-25
CATALOG #	G-300	G-302	G-301	G-303
ANTIMICROBIAL AGENT	20% EtOH			
STORAGE TEMPERATURE	2 - 8°C			

⁽¹⁾ Orientative values for binding capacity using BSA
 For laboratory use only. Not for use in diagnostic or therapeutic procedures.

Gold Biotechnology

St. Louis, MO

Web: www.goldbio.com

Ph: (314) 890-8778

email: contactgoldbio86@goldbio.com