

**AFFINITY HIS-TAG PURIFICATION  
TECHNICAL SPECIFICATIONS  
Empty Chromatography Columns**

<b>COLUMN MATERIAL</b>	Polypropylene	
<b>FRIT MATERIAL</b>	Polyethylene	
<b>FRIT PORE SIZE</b>	20µm	
<b>CHEMICAL STABILITY</b>	Stable in all commonly used reagents	
<b>CATALOG NUMBER</b>	<a href="#">P-301</a>	<a href="#">P-302</a>
<b>CAPACITY</b>	12 mL	35 mL

**For laboratory use only. Not for use in diagnostic or therapeutic procedures.**